

## Drug Recall

**DATE OF RECALL:** July 17, 2025

**DRUG NAME:** Dexcom G6 and G7 Receivers

**RECALLING FIRM:** Dexcom

**REASON FOR RECALL:** Dexcom announced a recall of its Dexcom G6 and G7 glucose monitoring receivers because a problem with the speaker may cause it to fail to make an alert sound when blood sugar is dangerously low or high.

**FDA LINK:** N/A

**OTHER DETAILS:** Do not use the recalled Dexcom receiver. Members with an affected receiver should contact Dexcom Technical Support at 1-844-478-1600 to request a free replacement.

**RECALLED PRODUCT:**

Product Description	NDC Number	Unique Device Identifier	Lot Number
Dexcom G6 Receiver Device	8627-0091-11	STK-FM-001 (MT27408-1) STK-FR-001 (MT27408-1)	00386270003232 00386270003225
Dexcom G7 Receiver Device	8627-0078-01	STK-AT-011(MT26403-0) STK-AT-012(MT26403-0) STK-AT-013(MT26403-0) STK-GT-001(MT26403-1) STK-GT-008(MT26403-5) STK-GT-013(MT26403-1) STK-GT-019(MT26403-4) STK-GT-100(MT26403-2) STK-GT-109(MT26403-2) STK-GT-113(MT26403-3)	00386270002198 00386270002518 00386270002525 00386270004291 00386270004260 00386270002273 04595986293080 00386270002020 00386270002129 00386270003638